

Rac1

Catalog Number: 26005

Gene Symbol: RAC1

Description: Anti-Rac1 Mouse Monoclonal Antibody

Background: Small GTPases are a super-family of cellular signaling regulators. Rac belongs to the Rho sub-family of GTPases that regulate cell motility, cell division, and gene transcription. In mammals, there are three Rac genes: Rac1, Rac2 and Rac3. GTP binding increases the activity of Rac, and the hydrolysis of GTP to GDP renders it inactive. GTP hydrolysis is aided by GTPase activating proteins (GAPs), while exchange of GDP for GTP is facilitated by guanine nucleotide exchange factors (GEFs).

Immunogen: Recombinant full length protein of Rac1

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA 1:1000-1:5000

WB 1:1000-1:2000

Concentration: 1 mg/ml

Host: mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid

Storage buffer:


Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rac1 antibody recognizes Rac1 from vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



Western blot analysis of Rac1 proteins. MEF whole-cell lysate was used for western blots. The monoclonal anti-Rac1 antibody was pre-incubated without (lane 1) or with (lane 2) the antigen protein (Rac1) to show the specificity.

For research use only. Not for diagnostic or therapeutic applications.