

G α_s

Catalog Number: 26006

Gene Symbol: Gnas

Description: G α_s Mouse Monoclonal Antibody

Background: Heterotrimeric G proteins are critical cellular signal transducers. G α_s represents one sub-family of G proteins that could mediate the stimulation of adenylyl cyclases.

Immunogen: Recombinant full length protein of G α_s

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA 1:1000-1:5000,

WB 1:1000-1:2000,

Host: mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid

Storage buffer:


Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: G α_s antibody recognizes G α_s of vertebrates

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



For research use only. Not for diagnostic or therapeutic applications.